

LIFEMED

T6100 / T6120

'S' Series Change

This single section couch is ideal for use in the care of disabled or dependent people and comes in a variety of lengths to suit your requirements. Supplied with two side rails for added safety (see optional accessories).

Features

- Two different choices of mid section length either 150cm or 180cm long
- Retractable rubber castors
- Adjustable foot for uneven floors
- Fire retardant vinyl upholstery



German Gas Struts & Motors




Specifications	150cm	180cm
Minimum Height	50cm	50cm
Maximum Height	96cm	96cm
Cushion Width	70cm	70cm
Head Section Length	N/A	N/A
Mid Section Length	150cm	180cm
Leg Section Length	N/A	N/A
Total Length	150cm	180cm
Lift Mechanism	Electric	Electric
Breathing Hole	0	0
Retractable Wheels	Standard x 4	Standard x 4
User Weight Rating	205kg	205kg


Warranty

Frame	3 Years
Electronic Components	3 Years
Upholstery	3 Years

Dimensions for 150cm and 180cm tables:



Height Range for 150cm and 180cm tables:



Code	Description	Optional Accessories:
MMC-METCHA150	T6100W "S" Series Change Table - length 150cm	STOCK
MMC-METCHA180	T6120W "S" Series Change Table - length 180cm	MMC-METPAPERHOLD Paper Towel Holder & Towel Roll
		NON-STOCK
		MMC-T6009M Chrome Side Rail
		MMC-T6005 Padded Side Rail