

LIFEMED T6718 'S' Series Massage

The Massage Table has a specially designed head piece that adjusts to suit the face of each patient. The head piece also tilts and the arm rests are located underneath and forward of the head cushion for patient comfort. The table facilitates excellent access to the patient's thoracic, and lower spinal regions.

Features

- The head rest has gas struts to help facilitate movement
- The wings provide a wider surface area beside the head rest for added comfort
- The footrest gives an extra 20cm when added, which is ideal for taller patients



German
Gas Struts
& Motors

*The wings and footrest come standard in Australia/New Zealand

Specifications

| | |
|---------------------|--------------|
| Minimum Height | 49cm |
| Maximum Height | 93cm |
| Cushion Width | 67cm |
| Head Section Length | 31cm |
| Mid Section Length | N/A |
| Leg Section Length | 150cm |
| Total Length | 184cm |
| Head Tilt Angle | -45°/10° |
| Lift Mechanism | Electric |
| Breathing Hole | 1 |
| Retractable Wheels | Standard x 4 |
| User Weight Rating | 205kg |

Warranty

| | |
|-----------------------|---------|
| Frame | 3 Years |
| Electronic Components | 3 Years |
| Upholstery | 3 Years |

Dimensions:






Table Positions:



Height Range:



| Code | Description | Optional Accessories: |
|----------------|--|--|
| MMC-METMESSAGE | T6718 'S' Series Massage Table - 67cm wide | NON-STOCK MMC-CU2115 Upholstery Head |