

LIFEMED

T8575 Elite Traction

The Elite Traction Table is a state-of-the-art traction table with a four piece top offering two friction free lockable sections as well as a lift up back section with a removable extension bar allowing for seated cervical traction. The design of this table allows for complete freedom in all traction set ups.

This table comes standard with an "Elevating Surround Bar", as shown in the picture.



German Gas Struts & Motors

Metron Traction Package



Specifications

Minimum Height	55cm
Maximum Height	101cm
Cushion Width	60cm
Head Section Length	41cm
Mid Section Length	41cm & 26cm
Leg Section Length	91cm
Total Length	202cm
Head Tilt Angle	-30°/30°
Leg Tilt Angle	0°/90°
Lift Mechanism	Electric
Breathing Hole	1
Retractable Wheels	Standard x 2
User Weight Rating	205kg

Warranty

Frame	3 Years
Electronic Components	3 Years
Upholstery	3 Years

Dimensions:





Table Positions:



Height Range:



Code	Description	Optional Accessories:
MMC-METETRAC	T8575 Elite Traction Table - 60cm wide	STOCK MMC-METPAPERHOLD Paper Towel Holder & Towel Roll MMC-METFLEXSTOOL Flexion Stool