

# LIFEMED

## T8714 Elite 3-Section

The Elite 3 is our most popular table with the 70cm becoming the standard in physiotherapy clinics. Featuring single gas strut configuration at head and foot section, "twin lift" stability, open base design, four levelling feet and available in a choice of colours, the Elite range combines quality and durability with style and elegance.

Complete and unrestricted access for a range of patient treatments.

Standard with face hole and plug. Comes with four wheels, as standard. Elevating surround bar available as an optional accessory.




### Specifications

Minimum Height	49cm
Maximum Height	96cm
Cushion Width	60cm, 70cm
Head Section Length	65cm
Mid Section Length	66cm
Leg Section Length	65cm
Total Length	202cm
Head Tilt Angle	-40°/80°
Leg Tilt Angle	-20°/80°
Lift Mechanism	Electric
Breathing Hole	1
Retractable Wheels	Standard x 4
User Weight Rating	249kg
Inclusions	Single Footswitch


### Warranty

Frame	3 Years
Electronic Components	3 Years
Upholstery	3 Years


### Dimensions:



### Table Positions:



### Height Range:



### Code

### Description

MMC-METELITE3	T8714	Elite 3-Section - 60cm wide
MMC-METELITE370	T8714W	Elite 3-Section - 70cm wide

### Optional Accessories

<b>STOCK</b>	
MMC-METPAPERHOLD	Paper Towel Holder & Towel Roll
MMC-T8008	Elevating Surround Bar
<b>NON-STOCK</b>	
MMC-T6009	Chrome Side Rail (for 60cm wide)
MMC-T6009X	Chrome Side Rail (for 70cm wide)