

LIFEMED

T8722 Elite 2-Section with Short Head



Designed and constructed using modern engineering techniques and manufactured from the finest materials to the highest standards.

Featuring single gas strut configuration at head section, "twin lift" stability, open base design and four levelling feet, the Elite range combines quality and durability with style and elegance.

In addition all models give complete and unrestricted access for a range of patient treatments.

Standard with face hole and plug. Comes with four wheels, as standard. Elevating surround bar available as an optional accessory.



Specifications

Minimum Height	50cm
Maximum Height	97cm
Cushion Width	60cm
Head Section Length	40cm
Leg Section Length	146cm
Total Length	200cm
Head Tilt Angle	-75°/15°
Lift Mechanism	Electric
Breathing Hole	1
Retractable Wheels	Standard x 4
User Weight Rating	249kg

Warranty

Frame	3 Years
Electronic Components	3 Years
Upholstery	3 Years

Dimensions:






Table Positions:



Height Range:



Code	Description	Optional Accessories:
MMC-METELITE2S	T8722 Elite 2-Section with Short Head - 60cm wide	STOCK MMC-METPAPERHOLD Paper Towel Holder & Towel Roll MMC-T8008 Elevating Surround Bar