

LIFEMED

T8725 Examination Chair MK2

Designed and constructed using modern engineering techniques and manufactured from the finest materials to the highest standards. Features "twin lift" stability, open base design, and four levelling feet. The MK2 offers three electric actuators, controlled by a foot switch. This controls the height, back rest and tilt angle adjustment. The leg rest has a smooth gas-strut adjustment. This allows perfect positioning for both patient and practitioner. In addition all models give complete and unrestricted access for a range of patient treatments.

Comes with four wheels, as standard.



German Gas Struts & Motors



Specifications

Minimum Height	59cm
Maximum Height	103cm
Cushion Width	55cm
Head Section Length	86cm
Mid Section Length	44cm
Leg Section Length	40/11cm
Total Length	202cm
Head Tilt Angle	0°/80°
Mid Tilt Angle	15°
Leg Tilt Angle	-40°/10°
Lift Mechanism	Electric
Breathing Hole	0
Retractable Wheels	Standard x 4
User Weight Rating	205kg

Warranty

Frame	3 Years
Electronic Components	3 Years
Upholstery	3 Years

Dimensions:





Table Positions:



Height Range:



Code	Description
MMC-METPOD2	T8725 Examination Chair MK2 - 55cm wide