

BEAUTY, SPA AND MASSAGE TABLE

A electric hi lo table is perfect for use in Beauty, Day Spa and Massage settings. It can be used in numerous comfortable positions with the head and leg section adjusted by gas struts.



FEATURES

- T motion motor for electric Hi Lo
- Lift height variation through footswitch
- Height variation from 57-98cm
- Adjustable backrest through gas operated strut to +85°
- Adjustable footrest through gas operated strut from +5 to -75°
- Adjustable feet and rounded legs for safety
- Leg extension
- Lift capacity of 225kg from the lowest point
- Retractable wheel system
- Luxurious padding for comfort
- Adjustable backrest
- Adjustable headrest
- Removable armrests

Available in white

Certification: CE, FDA/93, EEC, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz
sales@lifemed.co.nz
0800 893 952

ORDERING INFORMATION

Part Number	Size	Colour
BTW	186x50	White

DIMENSIONS

- Back section: 58cm
- Middle section: 45cm
- Leg section: 41cm
- Overall length: 186cm
- Adjustable Height: 48cm to 90cm
- Loading weight: 225Kgs

