

# METROCARE CHIROPRACTIC TABLE

Fully adjustable chiropractic couch with electric Hi Lo operation.



metrocare  
by LifeMed

## FEATURES

- T motion motor for electric Hi Lo
- Lift height variation through footswitch
- Minimum height of 54cms to maximum height of 87cm in 30 seconds
- Adjustable long board to +15°
- Height adjustable headrest from +25° to -20°
- 2 fixed wheels and 2 adjustable legs
- Leg extension
- Lift capacity of 225kg from the lowest point
- 5cm diameter wheels
- Designed with 4 manual drops:
  - Cervical
  - Thoracic
  - Lumbar
  - Pelvic

Available in navy blue and black

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

## ORDERING

Please contact us for all orders and for delivery options.

[www.lifemed.co.nz](http://www.lifemed.co.nz)  
[sales@lifemed.co.nz](mailto:sales@lifemed.co.nz)  
0800 893 952

## ORDERING INFORMATION

| Part Number | Size       | Colour    |
|-------------|------------|-----------|
| PM-CTNB     | 180-193x55 | Navy Blue |
| PM-CTBL     | 180-193x55 | Black     |

## DIMENSIONS

- Foot section length 13cm
- Long section length: 69.5cm
- Tapered section Length: 20.5 & 22.5cm
- Head section length: 40cm
- Upholstery width: 55cm
- Adjustable Height: 54cm to 87cm
- Loading weight: 225Kgs

