

CONTOURED HI LO ELECTRIC MASSAGE COUCH



Fully adjustable Metrocare contoured massage couch with an arm rest and cradle and electric Hi Lo operation by foot switch.



FEATURES

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section featuring double gas strut, angle adjustable from +45° to -45°
- Arm rest and cradle
- Adjustable foot for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable locking castors
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m2

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz
sales@lifemed.co.nz
0800 893 952

ORDERING INFORMATION

Part Number	Size	Colour
PM-CMNB	195x68	Navy Blue

DIMENSIONS

- Total Length: 197cm
- Body length: 165cm
- Head rest length: 32cm
- Upholstery width: 70cm
- Height: adjustable from 50cm to 90cm
- Open: 195 x 70 x (50-90)cm
- Loading weight: 250Kgs

