

# MULTIPURPOSE PODIATRY CHAIR

Metrocare 3 section podiatry couch with an electric Hi Lo operation by foot switch.



metrocare  
by LifeMed

## FEATURES

- Electric Hi Lo operation by foot switches
- Sleek design with rounded frame legs for safety
- Minimum height 56.5cm to maximum height 98cm in 30 seconds
- Dual motors. 1. Electric Hi Lo. 2. Electric mid-section
- Adjustable backrest featuring double gas struts, angle adjustable from 0° to +80°
- Adjustable seat section featuring one motor, angle adjustable from 0° to +25°
- Adjustable foot section featuring single gas strut
- Double removable arm supports
- Debris tray
- Able to lift up to 225kgs working weight from its lowest height
- Retractable castors 50mm for excellent mobility
- Table goes horizontal for people recovery in emergency situations
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m<sup>2</sup>

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

## ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz  
sales@lifemed.co.nz  
0800 893 952

## ORDERING INFORMATION

Part Number	Size	Colour
PM-PCNB	220x84	Navy Blue
PM-PCBL	220x84	Black

## DIMENSIONS

- Backrest length: 91cm
- Seat section length: 56cm
- Leg section length: 38cm
- Foot section length: 21cm
- Width: 84cms, 91cms inc arms
- Height: adjustable from 56.5cm to 98cm
- Open/flat: 220cm
- Loading weight: 225Kgs

