

# 2 SECTION CARDIOLOGY COUCH

Metrocare two section cardiology couch with comfortable sewn edge upholstery, central lock castors, electric backrest & two special cardiac cut outs allowing increased patient access.



## FEATURES

- Electric Hi Lo operation by foot switch (standard)
- Electric backrest for effortless operation
- Sleek design with rounded frame legs for safety
- Minimum height 47cm to maximum height 92cm in 30 seconds
- Adjustable head section featuring double gas strut, angle adjustable from flat to +75°
- Dual cut outs for easy patient access
- Cut outs operate manually with mechanical lock to hold in place
- 10 cm central locking castors for safety and mobility
- T Motion motor
- Comfortable sewn edge upholstery
- Padding thickness 5cm
- Foam Density: 36g/m2

Certification: CE, FDA, ISO9001-2008, ISO13485-2016

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

## ORDERING

Please contact us for all orders and for delivery options.

[www.lifemed.co.nz](http://www.lifemed.co.nz)  
[sales@lifemed.co.nz](mailto:sales@lifemed.co.nz)  
0800 893 952

## ORDERING INFORMATION

Part Number	Size	Colour
CARNB	195x66	Navy Blue

## DIMENSIONS

- Total Length: 195cm
- Width: 66cm
- Head section length: 67cm
- Leg section length: 123cm
- Height: adjustable from 58cm to 102cm
- Loading weight: 225Kgs

