

NEUROLOGICAL COUCH

Metrocare Designed for rehabilitation this strong two section neurological bobath has been designed with the therapist in mind. The large table top provides plenty of space for patient and therapist and three nose holes and plugs provides flexibility of treatment.



metrocare
by LifeMed

FEATURES

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 49cm to maximum height 93.5cm in 30 seconds
- Adjustable backrest section featuring double gas strut, angle adjustable to +65°
- Head section has three face holes and plugs
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m²

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz
sales@lifemed.co.nz
0800 893 952

ORDERING INFORMATION

Part Number	Size	Colour
NT120NB	190x120	Navy Blue

DIMENSIONS

- Total Length: 190cm
- Width: 120cm
- Backrest section length: 65cm
- Long section length: 121cm
- Height: adjustable from 49cm to 93.5cm
- Loading weight: 225Kgs

