

SINGLE SECTION CHANGE COUCH

Metrocare single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.



metrocare
by LifeMed

FEATURES

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Minimum height 48cm to maximum height 93.5cm in 30 seconds
- Face hole and plug
- Padded side rails
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m²

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz
sales@lifemed.co.nz
0800 893 952

ORDERING INFORMATION

Part Number	Size	Colour
PM-ET1NB	180x66	Navy Blue

DIMENSIONS

- Total Length: 180cm
- Width: 66cm
- Height: adjustable from 48cm to 93.5cm
- Loading weight: 225Kgs

