

# PROCEDURE CHAIR

The Metrocare procedure chair is a versatile multi modality chair that utilises fully electric functionality along with a strong column lift. It is perfect for the rigours of modern day use and will suit a number of different settings.



## FEATURES

- Electric Hi Lo operation by hand control
- Can achieve the Trendelenburg position
- Strong column lift
- Contoured back rest
- Gas strut adjustable face rest
- Debris tray
- Small footprint frame with rounded frame legs for safety
- Minimum height 59cm to maximum height 89cm in 30 seconds
- Adjustable feet for uneven surfaces
- 10cm lockable castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Comfortable foam padding

Certification: CE, FDA, ISO9001-2008, ISO13485-2016  
7 year warranty on frame and 3 year warranty on motor.

## ORDERING

Please contact us for all orders and for delivery options.

[www.lifemed.co.nz](http://www.lifemed.co.nz)  
[sales@lifemed.co.nz](mailto:sales@lifemed.co.nz)  
0800 893 952

## ORDERING INFORMATION

Part Number	Size	Colour
PRCHNB	198x64	Navy Blue

## DIMENSIONS

- Total Length: 198cm
- Width: 64cm
- Height: adjustable from 59cm to 89cm
- Adjustable back section: 0° to +85°
- Adjustable foot section: 0° to -85°
- Seat tilt: -5° to +35°
- Loading weight: 250Kgs

