

# VETERINARY LIFT TABLE

Metrocare height adjustable 1 section table with electric Hi Lo operation by foot bar (standard) for ease of operation.

This table is battery operated, eliminating the need for cables on the floor.

This single section table is specifically designed for veterinary use and features a simple design with a practical bench top in stainless steel.



## FEATURES

- Electric Hi Lo operation by foot bar (standard) or hand control (special order)
- Sleek design with rounded frame legs for safety
- Minimum height 41cm to maximum height 96cm in 30 seconds
- Adjustable feet for uneven surfaces
- 10cm lockable castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Durable stainless steel top
- Hygienic and practical for animal use
- Rechargeable batteries

Certification: CE, FDA, ISO9001-2008, ISO13485-2016  
7 year warranty on frame and 3 year warranty on motor.

## ORDERING

Please contact us for all orders and for delivery options.

[www.lifemed.co.nz](http://www.lifemed.co.nz)  
[sales@lifemed.co.nz](mailto:sales@lifemed.co.nz)  
0800 893 952

## ORDERING INFORMATION

Part Number	Size	Colour
ETVET	110x56	Stainless steel

## DIMENSIONS

- Total Length: 110cm
- Width: 56cm
- Height: adjustable from 41cm to 96cm
- Loading weight: 250Kgs

