

3 SECTION SHORT HEAD ELECTRIC COUCH

Metrocare fully adjustable 3 section couch with electric Hi Lo operation by foot switch.



FEATURES

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 47.5cm to maximum height 93cm in 30 seconds
- Adjustable head section featuring double gas strut, angle adjustable from -45° to +40°
- Adjustable leg section featuring double gas struts, angle adjustable to +80°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m²

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is stain, mildew, oil, water and wear resistant.

5 year warranty on frame and 3 year warranty on motor.

ORDERING

Please contact us for all orders and for delivery options.

www.lifemed.co.nz
sales@lifemed.co.nz
0800 893 952

ORDERING INFORMATION

Part Number	Size	Colour
PM-ET3NB	195x60	Navy Blue
PM-ET3BL	195x60	Black

DIMENSIONS

- Total Length: 195cm
- Width: 60cm
- Head section length: 40cm
- Middle section length: 54cm
- Leg section length: 93cm
- Height: adjustable from 47.5cm to 93cm
- Loading weight: 225Kgs

